


Disposable Set for  
**Percutaneous Trigeminal  
Nerve Compression**



Composition:

- > 01 Balloon of 1,5 mm;
- > 01 Interchangeable tips puncture;
- > 01 Button for puncture;
- > 01 Syringe of 1 ml;
- > 01 Cannula set guide;
- > 01 Mandrel with conical tip;
- > 01 Mandrel with angular tip.

**Technical Description:**

**Puncture:**

Material: Stainless steel  
Length - 220 mm  
Diameter - 1,7 mm  
Tips Cut - 15°

**Mandrel of the angular tip:**

Material: Stainless steel  
Length - 123,5 mm  
Diameter - 1,7 mm  
Tip Cut - 15°

**Balloon:**

Length - 400 mm  
Diameter - 1,5 mm  
Maximum capacity of liquid - 0,70 ml  
Maximum Strength of Traction (Insuf Balloon) - 1,5 lbs  
Diameter of the insufflate balloon - 9 mm

**Mandrel of the conical tip:**

Material: Stainless steel  
Length - 123,5 mm  
Diameter - 1,7 mm

**Syringe:**

Capacity - 1ml

**Cannula set guide:**

Material: Stainless steel  
Length - 100 mm  
Diameter - 2,1 mm

The images in this catalogue are illustrative.

KF0001FAA-REV.01-VER.ING 19/11/2009. All brands cited are the property of Micromar.

 **MICROMAR**

**CE**  
1023

**Micromar Europe S.L.**  
Avenida de Madrid nº19 - F  
ZIP 45003  
Toledo - Spain  
Phone/Fax +34 925 257 683  
europe@micromar.com  
www.micromar.com

**Micromar Ind. Com. Ltda.**  
Rua Orense nº396  
ZIP 09920-650  
Diadema - São Paulo - Brasil  
Phone/Fax +55 11 4057-5722  
sac@micromar.com  
www.micromar.com.br