


Pins for Cranial Fixation

Integration to the smallest details.

▲ Micromar Disposable Cranial Fixation Pins were especially designed with angled points, which guarantee the effectiveness of the procedure and assures less aggressive fixation for the patient. The pediatric pins also have a structure that prevents perforation beyond the safe fixation limit for children, providing more confidence to the professional. Micromar offers stainless-steel, aluminum, titanium and ceramic Cranial Fixation Pins, developed to meet different surgical needs. Also, all models are compatible with the main brands of cranial fixation systems on the market.



Color	Material	Code	Use
White	Ceramic	MD0130FAA	Adult
Green	Stainless-steel	MD0100FAA	Adult
Grey	Aluminum	MD0110FAA	Adult
Yellow	Titanium	MD0120FAA	Adult



Color	Material	Code	Use
White	Ceramic	MD0330FAA	Pediatric
Green	Stainless-steel	MD0300FAA	Pediatric
Grey	Aluminum	MD0310FAA	Pediatric
Yellow	Titanium	MD0320FAA	Pediatric

*Illustrative chart representing the interference level of each material on image acquisition.

The images of this catalogue are illustrative.
REV:02 VER:ENG 12/02/2010. All marks cited at this catalogue are Micromar property.

MICROMAR

CE
1023

Micromar Europe S.L.
Avenida de Madrid nº19 - F
ZIP 45003
Toledo · Spain
Phone/Fax +34 925 257 683
europe@micromar.com
www.micromar.com

Micromar Ind. Com. Ltda.
Rua Orense nº396
ZIP 09920-650
Diadema · São Paulo · Brasil
Phone/Fax +55 11 4057-5722
sac@micromar.com
www.micromar.com.br